Changing to a bigger 215g tin from February 2025

Same Formulation. Same Scoop Size.



More Scoops Per Tin Means Fewer Tins Per Prescription

Resource ThickenUp Clear tin size is updating on the GP prescribing systems from February 2025. The 127g tin will be discontinued and replaced with a 215g tin, which will have a new PIP code.

Old Case Amount	Old PIP code	127g Tin Price	NEW Case Amount	NEW PIP code	215g Tin Price
Resource® ThickenUp® Clear 12 x 127g	355-4433	£9.02	Resource® ThickenUp® Clear 6 x 215g	433-7291	£15.27





To request your **FREE sample** please scan the QR code. www.nhscsampleservice.co.uk

For more information visit www.nestlehealthscience.co.uk/brands/resource/215g-tin or call our customer care line on 00800 6887 4846



Resource® ThickenUp® Clear (215g)

Prescribing guide† for new 215g Tin.

[†]To be used as a general guide, for Resource *Thicken Up * Clear only. Based on the recommended usage of Resource *Thicken Up * Clear. **Same scoop size**, **Same product**, **New tin size**.

Resource ThickenUp Clear is a transparent and tasteless, instant food and drink thickener, designed to support the hydration of patients with dysphagia (swallowing difficulties).



1DAY	LEVEL 1 SLIGHTLY THICK		LEVEL 2 MILDLY THICK		LEVEL 3 MODERATELY THICK		LEVEL 4* EXTREMELY THICK	
ml	Scoops	g	Scoops	g	Scoops	g	Scoops	g
200	1	1.2	2	2.4	4	4.8	6	7.2
600	3	3.6	6	7.2	12	14.4	18	21.6
1000	5	6	10	12	20	24	30	36
2000	10	12	20	24	40	48	60	72

7 DAYS	LEVEL 1 SLIGHTLY THICK		LEVEL 2 MILDLY THICK		LEVEL 3 MODERATELY THICK		LEVEL 4* EXTREMELY THICK	
ml	NEW Tins [‡]	Scoops	NEW Tins [‡]	Scoops	NEW Tins [‡]	Scoops	NEW Tins [‡]	Scoops
200	1	7	1	14	1	28	1	42
600	1	21	1	42	1	84	1	126
1000	1	35	1	70	1	140	2	210
2000	1	70	1	140	2	280	3	420

28 DAYS	LEVEL 1 SLIGHTLY THICK	LEVEL 2 MILDLY THICK	LEVEL 3 MODERATELY THICK	LEVEL 4* EXTREMELY THICK	
ml	NEW Tins [‡]	NEW Tins [‡]	NEW Tins [‡]	NEW Tins [‡]	
200	1	1	1	1	
600	1	1	2	3	
1000	1	2	4	5	
2000	2	4	7	10	

1 scoop =1.2 g. For food, add Resource® ThickenUp® Clear as needed to achieve the appropriate consistency. *For level 4 extremely thick liquids, evaluated IDDSI compliant based on spoon tilt test. It is the responsibility of the person administering Resource® ThickenUp® Clear to ensure that the liquid or food is mixed to the appropriate consistency. ‡ New 215g Resource® ThickenUp® Clear tin.

The tables above help determine the amount of Resource® ThickenUp® Clear required to thicken drinks for a patient during a given period. Fluid intake and individual patients' consistency requirements should be determined by a healthcare professional. It is essential that recommended consistency requirements are followed.

